





Abdominal Pads

SS 102 A

-  Abdominal Pads
-  Pastillas Abdominales
-  Pads abdominais
-  compresses abdominales



- Used in surgical applications
- Made from neatly stitched, fine absorbant, non-woven gauze fabric, in multiple layers
- X-ray thread for ease in visualization during scanning
- Manufactured as per British Pharmacopia (BP)/USP guidelines
- Option of 4 ply, 8ply or 12ply

Packing

- Available in Sterile/Non Sterile pack.
- Individually packed or available in combi pack of 2, 4 or 100pcs.



Gross: 9kg or customized.





1000 pcs or customized

SIZES

25cm x 25cm, 30 cm x 30 cm & 35cm x 35cm

Surgical Pad/Gauze & Cotton Tissue /Combination Dressing

SS 1008

-  Surgical Pad/Gauze & Cotton Tissue/Combination Dressing
-  Pad Quirúrgica/Gasa Y Tejidos De Algodón/Vestidor Combinación
-  Pad cirúrgico / Gauze & Cotton Tissue / Molho Combinação
-  Pad chirurgicale / Gaze & Coton Tissu / Combinaison Dressing



Packing

- Sterile, individually packed in peelable pouch pack.
- Carton of 500pcs.



Gross: 7.32kg



280 x 482 x 788mm

SIZES

6cm x 6cm, 10cm x 10cm & 10cm x 20cm

- Used for dressing and general medical applications
- Made from fine absorbent cotton covered with a layer of non-woven gauze fabric
- Manufactured as per British Pharmacopia (BP)/USP guidelines