


## Silicon Foley Balloon Catheter (BH Model)

2 WAY | SMD 517  
3 WAY | SMD 519

- All Silicon Foley Balloon Catheter
- Todo Es Silicio Foley Globo Sonda
- Todos cateter de silicone Foley Balloon
- Toutes les cathéter Foley Silicon Balloon



High elastic tear resistant balloon



Flow Control Clamp

- Provided with flow control clamp to ensure proper use of bladder muscles .
- Self lubricated silicon shaft provides low friction surface to allow trauma free removal.
- Balloon made from low durometer, highly elastic silicon which provides very good tear strength.

TWO WAY		THREE WAY	
Size in FR	Balloon Capacity	Size in FR	Balloon Capacity
6	3cc	16,18, 20, 22, 24 & 26	30cc - 50cc
8 -10	3cc - 5cc		
12 - 14	5cc - 15cc		
16 - 26	30cc - 50cc		

### Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 10; Master Carton of 600



Gross: 15.3kg (14 FR)



690 x 330 x 530mm

SIZES										
FR	8	10	12	14	16	18	20	22	24	26
Colour Code	Black	Grey	White	Green	Orange	Red	Yellow	Voilet	Blue	Pink