





## Urine/Stool Specimen Container

SMD 902

-  Urine / Stool Specimen Container
-  Contenedor De Muestras De Orina / Heces
-  Urina / Stool Specimen Container
-  Urine / échantillon de selles Conteneurs







- Used for urine/stool sample collection
- Made from medical grade translucent polypropylene/clear polystyrene
- Available in sterile, individual tearpouch pack



## Urine Specimen Container

SMD 901

-  Urine Spesimen Container
-  Contenedores De Muestras De Orina
-  Urine Specimen Kontainer
-  Échantillon d'urine Conteneurs

### Packing

- Pack of 25. Carton of 600



Gross: 9.23kg








660 x 432 x 508mm

- Used for urine sample collection
- Graduated air-tight container available in 30/70ml capacity

## Blood Lancet

SMD 900

-  Darah Lancet
-  Lanceta De Sangre
-  Lancet sanguie
-  Lancette



- Used as a manual blade to make a limited skin incision for blood collection
- Made from medical grade stainless steel and PVC
- Ultrasonically cleaned, with super sharp tip for the instant incision

### Packing

- Sterile, individually packed in peelable pouch pack.
- Pack of 200. Carton of 20000pcs.



Gross: 13kg



432 x 432 x 432mm