

## Chest Drainage System /Under Water Sealed Drainage System

SMD 1043

- Dada Sistem Drainase / Under Water Sealed Sistem Drainase
- Pecho Sistema de Drenaje / Under Water Sistema de Drenaje Sellado
- Peito Sistema de Drenagem / Under Water sistema de drenagem estanque
- Système de drainage thoracique / Système de drainage scellés sous l'eau



- Used to remove air (Pneumothorax) or fluid (Pleural Effusion, blood) or pus (empyema) from the intrathoracic space
- Made from non-toxic medical grade PVC
- Bottle with graduation markings to determine precise drain volume
- Clear initial level markings for easy and correct usage
- Suction port to connect bottle with suction
- Universal connectors and kink-resistant tubings for smooth unrestricted flow
- Hangers available for ease of use
- Available in 500ml capacity (single chamber), 1000ml capacity (single chamber) and 2000ml (double chamber)

SIZES	
Single Chamber	500ml
Double Chamber	2000ml

### Packing

- Individually packed in Inner Carton.
- Master Carton of 18 pcs. (500ml)
- Master Carton of 10 pcs. (1000ml)
- Master Carton of 10 pcs. (2000ml)



Gross: 11kg (2000ml)



770 x 360 x 450mm