



Clearpore

Transparent Perforated Plastic Film Tape
SS 106 C


- Transparan Tape berlubang Plastik Film
- Cinta de película plástica transparente perforada
- Fita de filme plástico transparente perfurada
- Ruban de film plastique transparent perforé



Breathable



Water resistant



Bi-directional tearing

- Made from perforated plastic film for good breathability
- Uniformly coated with hypoallergenic adhesive
- Strong adhesion water resistant
- Painless removal
- Perfect bidirectional tearing

Packing

- Carton of 60 Boxes

SIZES & PACKING		
Size	Carton Size	Gross Weight
1/2", 1", 2" x 9.1mtr	480x430x285mm	15.0kg
3" x 9.1mtr	430x330x420mm	14.5kg
1/2", 1", 2" x 5.0mtr	390x360x285mm	9.0kg
3" x 5.0mtr	435x270x345mm	8.5kg