






Infant Feeding Tube Plain / Elastomer Coated

With Graduation | SMD 603 G
 Without Graduation | SMD 603
 Elastomer Coated | SMD 603 E



PHTHALATE FREE

-  Infant Feeding Tube
-  Tubo De Alimentación Infantil
-  Infant tubo de alimentação
-  Infant Feeding Tube



- Used for neonates and paediatric feeding
- Made from non-toxic, medical grade PVC
- Smooth coned distal end for atraumatic intubation, with two lateral eyes for efficient drainage
- Markings at 20cm and radio opaque line for accurate placement and X-ray visualization
- Proximal end with luer mount for connection with feeding device
- Marking from the tip 5 cm to 40 cm

Packing

- Sterile, Individually packed in peelable blister pack.
- Box of 200; Master Carton of 2000.

EFFECTIVE LENGTH	TOTAL LENGTH
50cm	52cm

SIZES							
FG	4	5	6	7	8	9	10
Colour Code	White	Grey	Green	White	Blue	Orange	Black

 Gross: 15.5kg (6 FG)  680 x 305 x 650mm