


## Kehr's "T" Tube

SMD 606

- Kehr's "T" Tube
- Tubo De Kehr "T"
- Tube "T" de Kehr
- "T" Tube de Kehr

- Used for short term post operative drainage of common bile duct
- Made from soft non-toxic, medical grade PVC
- Radio opaque line provided throughout the tube for X-ray visualization
- Super smooth tube surface to prevent clogging



PHTHALATE FREE

### Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 50; Master Carton of 500 (10 to 16 FG)
- Box of 50; Master Carton of 600 (18 to 20 FG)

Gross: 12.3kg (16 FG)

10 to 16 FG - 790 x 480 x 280mm  
18 FG - 680 x 480 x 520mm

SIZES	
Size	10FG , 12FG , 14FG , 16FG & 18FG
Length	20 x 60cm



PHTHALATE FREE



## Rectal Catheter

SMD 605

- Kateter dubur
- Sonda Rectal
- cateter rectall
- cathéter rectal

- Used for introduction of enema solution into rectum
- Made from non-toxic, non-irritant medical grade PVC
- Open distal end is open with one lateral eye for efficient drainage
- Proximal end with colour coded funnel connector for easy connection to the suction equipment
- Funnel connectors are colour coded for instant size identification
- Optional: Available in 40cm length

EFFECTIVE LENGTH	TOTAL LENGTH
37cm	40cm

### Packing

- Sterile, Individually packed in peelable pack.
- Box of 100; Master Carton of 1000.

Gross: 25.5kg (18 FG) 680 x 305 x 650mm

SIZES								
FG	12	14	16	18	20	22	24	26
Colour Code	White	Green	Orange	Red	Yellow	Voilet	Blue	Pink