

Nelaton / Urethral Catheter

SMD 502

- Nelaton / Urethral Cateter
- Nélaton / Urétrhal Catéter
- Nelaton / Urethral Catheter
- Nélaton / Urétrhal Cathétér



- Used for short term bladder catheterization
- Made from non-toxic, non-irritant, medical grade PVC
- Radio opaque line provided through out the catheter for X-ray visualization
- Distal end with lateral eyes and proximal end funnel connector, for maximum drainage and easy connection to the urine bag

Packing

- Sterile, Individually packed in peelable blister pack.
- Box of 100; Master Carton of 1000.

VERSION	EFFECTIVE LENGTH	TOTAL LENGTH
Male	370mm	400mm
Female	170mm	200mm

Gross: 14.2kg (14 FG) 680 x 305 x 650mm

SIZES								
FG	6	8	10	12	14	16	18	20
Colour Code	Light Green	Blue	Black	White	Green	Oragne	Red	Yellow



Male External Catheter Penile Sheath

SMD 503

- Male External Catheter Penile Sheath
- Catéter Externo Masculino Con Funda Pene
- Cateter masculino externo peniana Bainha
- Homme externe cathéter gaine pénienne

- Used for urine drainage in incontinent male patients
- Made from natural latex rubber
- Provided with adhesive strip for proper fixing of the catheter

Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 50; Master Carton of 500.

Gross: 7.1kg (Large) 640 x 410 x 220mm

SIZES	
Small	20mm
Medium	25mm
Large	30mm
Extra Large	35mm