





Nebulizer Mask

Child/Adult
SMD 710

-  Nebulizer Masker
-  Máscara Nebulizadora
-  Μάσκαρα Νεβυλίζατορ
-  Masque De Nébulisation



- Specially designed nebulizer chamber with nebulization rate of approximately 3 cc in 10 minutes.
- Average particle size generation is approximately 2 microns.
- Complete kit consists of: Nebulizer mask, Nebulizer chamber and 210 mm Star Lumen Tube.
- Sizes : Adult and paediatric

Packing

- Sterile, individually packed in blister pack
- Box of 10, Master Carton of 200.



Gross: 24.2kg



680 x 480 x 520mm