

## P.E. Pressure Monitoring Line

Spiral

SMD 715

- P.E. Pressure Monitoring Line Spiral
- P.E. Línea Espiral De Monitoreo De Presión
- P.E. Linha de Monitorização da Pressão
- P.E. Surveillance de la pression de ligne

- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from Kink resistant PE Tubing
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Internal Diameter 1.0 mm, Outer Diameter 2.0 mm

### Packing


- Sterile, Individually packed in peelable pouch & blister pack.
- Box of 100; Master Carton of 500



Gross: 13.5kg



790 x 480 x 280mm



| SIZES        |       |        |        |
|--------------|-------|--------|--------|
| Length       | 110cm | 150cm  | 200cm  |
| Prime Volume | 0.9ml | 1.40ml | 1.80ml |



## Pressure Monitoring Line

SMD 706

- Pressure Monitoring Line
- Línea De Monitoreo De La Presión
- Linha de Monitorização da Pressão
- Surveillance de la pression de ligne

- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from non-toxic medical grade PVC
- Can withstand a pressure of 580 PSI
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Optional lengths available

### Packing

- Sterile, Individually packed in peelable pouch & blister pack.
- Box of 50; Master Carton of 500



Gross: 14.5kg (200cm)



670 x 300 x 480mm

| SIZES        |       |        |        |
|--------------|-------|--------|--------|
| Length       | 110cm | 150cm  | 200cm  |
| Prime Volume | 0.9ml | 1.40ml | 1.80ml |