




Red Rubber Nelaton Catheter / Jaques IR Catheter

SMD 522

-  Red Karet Nelaton Kateter
-  Catéter de Nelaton de caucho rojo
-  Cateter vermelho de borracha Nelaton
-  Cathéter de Nelaton en caoutchouc rouge



- Red natural latex rubber catheter used across the spectrum of health care for all patient types
- High quality construction ensure consistency of use and patient comfort
- Smooth tip and funnel end made from a radiopaque material
- Two large opposing eyes for optimum drainage
- Suitable for male or female catheterization
- Product type - Intermittent Catheter
- Product style - Straight tip

Packing

- Sterile, Individually packed in peelable blister pack.
- Box of 10; Master Carton of 600.



Gross: 15kg



680 x 330 x 530mm

SIZES								
FR	8	10	12	14	16	18	19	20