





Steri-ADK

Abdominal Drainage Kit
SMD 1011

-  Perut Drainase Kit
-  Kit Drenaje Abdominal
-  Kit de Drenagem abdominal
-  Kit de drainage abdominale

- Used for abdominal drainage
- Made from non-toxic, medical grade PVC
- Smooth and round distal end with large atraumatic eyes for efficient drainage
- Radio-opaque line provided on the catheter for x-ray visualization
- Marked at every 2cm from the cast eye to ascertain the depth of placement
- Proximal end provided with universal taper connector for easy connection to drainage bag
- 2000ml capacity bag with graduated scale

Packing

- Box of 10; Master Carton of 100.



Gross: 9.7kg (24 FG)







720 x 340 x 360mm



Steri-Drain

Chest Bag / Water Sealed Drainage System
SMD 1012

-  Chest Bag / Water Sealed Drainage System
-  Bolsa De Pecho/ Sistema De Sellado Pra El Agua Del Drenaje
-  Saco Peito / Água Selado Sistema de Drenagem
-  Poitrine Sac / Eau Sealed système de drainage

Packing

- Individually packed in a polybag.
- Box of 10; Master Carton of 100pcs.

- Used for collection of wound and empyema
- Made from non-toxic, non-irritant PVC
- 1000ml capacity bag with graduated scale



Gross: 11kg



540 x 430 x 350mm