



## Steriderm

I.V Transparent

SS 116

-  Steriderm I.V Transparan
-  Steriderm I.V Transparente
-  Steriderm i.v Transparente
-  Steriderm i.v Transparent


- Used for Fixation of I.V Cannula / Catheter
- Made from Thin Polyurethane film and available in two variants -Transparent and Non-woven
- Breathable and water proof, with reinforced notches that reduce mechanical stress from heavy catheters
- Thin film ensures easy fixation of devices with contoured surfaces



### Packing

- Inner Box of 50 pcs.; Master Carton of 5000 pcs.





 Gross: 11kg

 → 437 x 360 x 225mm

## Steriderm I.V

Elastic Adhesive Base

SS 116 A

-  Steriderm I.V Elastis Adhesive Dasar
-  Steriderm I.V Adhesivo Con Base Elastica
-  Steriderm i.v Elastic Base Adesiva
-  Steriderm i.v Elastic base adhesive




- Used for Fixation of I.V Cannula / Catheter
- Made from high-strength and extra breathable elastic fabric
- Uniformly coated with hypoallergenic, anti-fungal, non-irritant adhesive

### Packing

- Sterile packed in peelable pouch pack.
- Inner box of paper - 100 pcs. Outer Carton of 4000 pcs.
- Inner box of plastic - 150 pcs. Outer Carton of 1800 pcs.

 Gross: 11.1kg

 → 635 x 483 x 483mm (4000 pcs.)  
570 x 490 x 290mm (1800 pcs.)