


Steriflon

I.V Cannula With Injection Port
SMD 1100

- Steriflon I.V Cannula With Injection Port
- Cánula De Steriflon I.V. Con Puerto De Inyección
- I.v Cannula com injeção Porto
- I.v canule à injection Port



- Used for intra venous infusion of medicinal fluids
- Made from non-toxic medical grade PTFE
- Formulated cannula tip and smooth surface for atraumatic catheterization
- Super sharp needle with two cutting edges
- Provided with injection port with non-return valve for extra medication

Packing

- Sterile, individually packed in peelable blister pack.
- Box of 100 pcs.; Master Carton of 1000 pcs.

Gross: 6.7kg

698 x 320 x 178mm

SIZES						
Size in FG	14	16	18	20	22	24
Colour Code	Orange	Grey	Green	Pink	Blue	Yellow
Length	45mm	45mm	45mm	32mm	25mm	19mm