

## Thoracic Trocar Catheter

SMD 1001

- Thoracic Trocar Kateter
- Catéter Trocar Torácica
- Trocar torácica Cateter
- Thoracique trocarts catheter

- Used for rapid non-operative plueral and chest drainage
- Made from non-toxic, medical grade PVC
- Catheter is snugly fixed over the trocar, and provided with large smooth eyes for efficient drainage
- Smooth taper shape tip
- Markings at 5, 10, 15 and 20cm from the distal end to measure depth of placement
- Proximal end fitted with taper shape connector for easy connection to drainage bottle/bag



PHTHALATE FREE



### Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 10; Master Carton of 100.

Gross: 8.5kg (28 FG)

550 x 400 x 270mm

SIZES										
Size	8	10	12	14	16	20	24	28	32	36
Length (mm)	230	230	230	230	230	390	390	390	390	390



PHTHALATE FREE



## Thoracic / Chest Drainage Catheter

SMD 1002

- Thoracic / Chest Drainage Catheter
- Catéter De Drenaje Torácico
- Thoracic / Chest drenagem por cateter
- Thoracique / Chest cathéter de drainage

- Used for post-operative drainage after cardio thoracic surgery
- Made from non-toxic, medical grade PVC
- Smooth and round open distal end with large eyes for maximum drainage
- Catheter is marked at every 2 cm from the last eye to ascertain the depth of placement
- Proximal end fitted with taper tongue connector for better forceps grip and smooth penetration
- Optional: Also available in right angle

### Packing

- Sterile, Individually packed in peelable pouch.
- Box of 20; Master Carton of 200

Gross: 11.5kg (28 FG)

680 x 305 x 650mm

SIZES									
Size	12	14	16	20	24	28	32	36	
Length (mm)	450	450	450	450	450	450	450	450	450