

## Tracheostomy Tube - Cuffed & Uncuffed

Cuffed | SMD 702 C  
 Uncuffed | SMD 702 P

- Trakeostomi Tube - diborgol & Uncuffed
- La cánula de traqueostomía - Abofeteado y sin manguito
- Traqueostomia - Cuffed & uncuffed
- Trachéotomie - revers et sans ballonnet



PHTHALATE FREE



Lockable obturator with provision for guidewire introduction

- Used to maintain a patent airway through tracheostomy and mechanically ventilate patients
- Made from non-toxic, non-irritant, medical grade clear PVC
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- **Lockable obturator** with facility to insert guide wires to facilitate intubation using the widely accepted Seldinger guide wire technique
- Anatomically designed tube and obturator with extra smooth and curved distal end for atraumatic intubation and extubation
- Smooth inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use

### Packing

- Sterile, Individually blister packed.
- Master carton of 50.



Gross: 4.75kg



580 x 240 x 290mm

### SIZES

3.5 mm, 4.0 mm, 4.5 mm, 5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm, 9.0 mm