

Umbilical Cord Clamp

SMD 801

- Umbilical Cord Clamp
- Abrazadera Cordón Umbilical
- Cordão Umbilical Braçadeira
- Pince cordon ombilic



- Used for clamping the umbilical cord of newly born babies
- Made from non-toxic, medical grade PVC
- Permanent locking system to prevent clamp reopening
- Finger grip to facilitate easy holding and locking of the clamp

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 100; Master Carton of 2000.



Gross: 10.5kg




670 x 300 x 480mm

Umbilical Catheter

SMD 1103

- Kateter Umbilical
- Catéter Umbilical
- cateter umbilical
- cathéter ombilical



- Used to access the umbilical vessels of pre-mature or newly - born babies
- Made from non-toxic, medical grade PVC
- Smooth round tip for atraumatic insertion and well-finished surface facilitate smooth passage in the vein
- Markings at every cm to ascertain depth of insertion, radio opaque line ascertain position during scans
- Proximal end provided with Luer female connector for easy connection to the I.V. therapy equipment

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 100; Master Carton of 1000



Gross: 8.5kg



680 x 305 x 650mm

SIZES

Total Length - 40cm

SIZES

FG	4	5	6	8	10
Colour Code	White	Grey	Green	Blue	Black