

Urine Collecting Bag With Measured Volume Meter

SMD 510

- Urine Collecting Bag With Measured Volume Meter
- Bolsa Colectora De Orina Con Medidor De Volumen
- Urine Bag Coletando Medido com medidor de volume
- Collecteur d'urine Sac Avec Mesuré indicateur de volume



- Used for measurement of hourly urine output
- 2000ml capacity bag made from non-toxic, non-irritant, medical grade PVC
- Rigid 250ml capacity bottle for accurate measurement of urine output
- Bottle directly connected to the bag to avoid contamination
- Push/Pull type bottom outlet for convenient emptying of urine bag
- Kink resistant 90cm tube with universal connector

Packing

- Sterile, Individually packed in peelable blister pack.
- Box of 1; Master Carton of 50.

Gross: 5.5kg 680 x 480 x 520mm

SIZES	
Bottle Capacity	250ml (Adult)
Bag Capacity	2000ml
Tube Length	90cm



Urine Collecting Bag With Hanger

SMD 504

- Urine Mengumpulkan Bag Hanger Dengan
- La orina que recoge el bolso con la suspensión
- Coleta de urina Saco Com Hanger
- Collecteur d'urine Sac Avec Hanger

- Used for short and long term urine collection drained through a catheter
- 2000ml capacity bag made from medical grade soft PVC
- Non-return valve to avoid back flow of Urine
- Kink resistant 90cm tube with universal connector for easy connection to the drainage catheter
- Graduated scale to measure urine output accurately
- Optional - top bottom outlet and molded handle for easy handling and carrying
- Available without hanger (SMD 504T)

Packing

- Sterile, Individually packed in Tear H.M pack.
- Box of 25; Master Carton of 400.

Gross: 25.2kg 680 x 480 x 520mm

SIZES	
Bottle Capacity	200ml
Tube Length	90cm