





Sterimed Surgical Gloves

SMD 1000

-  Sarung tangan bedah
-  Guantes Quirurgicos
-  luvas para cirurgia
-  Gants chirurgicaux



- Used for protection during general medical/surgical procedures
- Made from natural latex rubber
- Firm beaded cuffs, micro rough finish for better grip
- Bent fingers for greater comfort
- Lightly powdered

Packing

- Carton of 1000 pairs. Sterile



Gross: 28.3kg



483 x 660 x 457mm

SIZES

6.0, 6.5, 7.0, 7.5, 8, 8.5