





URIGATE™ - CATHETER VALVE

SMD 529

-  Urigate - Kateter Katup
-  Urigate - Válvula de catéter
-  Urigate - Válvula Cateter
-  Urigate - Valve de cathéter



Urigate inserted in Foley Catheter

- Suitable for use with both urethral and suprapubic catheters.
- Helps maintain bladder tone and capacity by promoting regular flushing.
- Designed to minimize the risk of blockage for uninterrupted use.
- Features an easy-to-operate mechanism, allowing for simple bladder emptying at the user's convenience.
- Enables patients to manage their own bladder function, supporting autonomy and mobility.
- Built for durability and safety, ensuring consistent performance and peace of mind.

Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 20; Master Carton of 400



Gross: 7.00kg



670 x 300 x 480mm