

Sterimed Preformed Endotracheal Tube - Cuffed & Plain

Conforming to ISO 5361

South Oral & North Nasal

North Nasal Cuffed | SMD 717 C

North Nasal Plain | SMD 717

South Oral Cuffed | SMD 718 SC

South Oral Plain | SMD 718 S



DEHP FREE

-  Preformed Endotracheal Tube - Cuffed & Plain
-  Tubo Endotraqueal Preformada - Con Puño/ Plano
-  Endotracheal pré-formada Tube - Cuffed & Plain
-  Endotrachéale préformes Tube - Revers & Plain



- Used for maxillofacial surgeries.
- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade PVC
- Preformed tubes reduce the risk of kinking and improve access to the operative field
- Tube can be temporarily straightened for ease of suctioning
- Smooth polished murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line
- High volume / low pressure cuff ensures efficient low pressure cuff seal (For cuffed tubes)
- Universal 15 mm connector at proximal end
- E.T.O sterile and pyrogen free. For single patient use

Packing

- Sterile, Individually peelable pouch packed.

SIZES - ORAL & NASAL												
ID (mm)	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5
O.D. (mm)	4.2	4.8	5.5	6.2	6.9	7.6	8.2	8.9	9.6	10.3	10.9	11.6

	Inner Box	Master Carton	Carton Dimensions	Gross Weight
South Oral Cuff/Plain	10 pcs.	200 pcs.	 660 x 320 x 380 mm	 7.56kg Approx.
North Nasal Cuff/Plain	10 pcs.	180 pcs.	 660 x 380 x 500 mm	 9.86kg Approx.