

FREEVENT®
Nasopharyngeal Airways
SMD 720

- Nasofaríngeo Airways
- Nasopharyngeal Airways
- nasofaríngeas Airways
- nasopharyngés Airways

- Used to secure and open airways for patients in case of emergencies
- Provided **with rotatable flange** for easy securing and pulling out.
- Made from soft, non-irritant medical grade PVC
- Anatomically designed with smooth round edges, thin wall and kink resistance for maximum patient comfort
- E.T.O. sterile & pyrogen free. For single use only

Packing

- Sterile individually peelable pouch/blister packed.
- Box of 25; Master Carton of 200 pcs.

Gross: 6.66kg 660 x 320 x 380mm



SIZES
3.0 mm, 3.5 mm, 4.0 mm, 4.5 mm, 5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm, 9.0 mm

Sterimed Guedel Airways / Oropharyngeal Airways
SMD 705

- Guedel Airways / Orofaríngeo Airways
- Cánulas De Guedel / Orofaríngeos Airways
- Guedel Airways / orofaríngeicos Airways
- Guedel Airways / oropharyngés Airways



- Used to maintain an unobstructed oropharyngeal airway
- Made from rigid, non-toxic medical grade LDPE
- Colour coded in-built bite blocks to avoid airway occlusion and biting of tongue

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 50; Master Carton of 400.

Gross: 8kg (Size 03) 000 to 5 - 770 x 360 x 450mm

SIZES			
Size	Colour Code	Size	Colour Code
000	-	2	Green
00	Blue	3	Yellow
0	Black	4	Red
1	White	5	Orange