



## Sterimed Nebulizer Mask

Child/Adult

SMD 710

-  Nebulizer Masker
-  Máscara Nebulizadora
-  Máscara Nebulizador
-  Masque De Nébulisation



DEHP FREE

- Specially designed nebulizer chamber with nebulization rate of approximately 3 cc in 10 minutes.
- Average particle size generation is approximately 2 microns.
- Complete kit consists of: Nebulizer mask, Nebulizer chamber and 210 mm Star Lumen Tube.
- Sizes : Adult and paediatric

### Packing

- Sterile, individually packed in poly pack
- Box of 20, Master Carton of 200.



Gross: 21.9kg



740 x 425 x 425mm

## Sterimed Nasal Oxygen Cannula

Twin Bore Oxygen Delivery Set

SMD 703 A

-  Nasal Oxygen Cannula Twin Bore Oxygen Delivery Set
-  Canula Nasal De Oxigeno Entrega Con Doble Juego
-  Nasal cânula de oxigênio Bore Twin Set entrega de oxigênio
-  Canule nasale d'oxygène Bore Twin Set oxygène de livraison



DEHP FREE

| SIZES  | TUBE LENGTH |
|--------|-------------|
| Adult  | 210cm (7ft) |
| Child  | 210cm (7ft) |
| Infant | 210cm (7ft) |

Total Length 250cm

\* (Length can be customized)

### Packing

- Non-Sterile, individually packed in a polybag.
- Box of 25, Master Carton of 250.



Gross: 16.6kg



790 x 480 x 280mm