

Sterimed Suction Catheters

Plain / Normal | SMD 700 P
 Thumb Control | SMD 700 ATC
 Finger Tip Control | SMD 700 BFT
 Tiemann Tip | SMD 700 T

 Suction Catheters
 Catéteres De Succión
 Cateteres de sucção
 Los catéteres de succión

CE
 0123



Tip & Eye Configurations



Tiemann Tip



Double Eye



Single Eye

- Used for removal of secretions from the patient's mouth, trachea and bronchial tube
- Made from non-toxic, non-irritant medical grade PVC
- Lateral eye at distal end for effective suctioning
- Proximal end of tube provided with universal connector for easy connection with suction source
- Optional: Thumb controlled, Finger tip controlled, Atraumatic tip
- Available with graduation

Packing

- Sterile, individually packed in blister pack
- Box of 100; Master Carton of 1000.



Gross: 15.2kg (14 FG)



680 x 305 x 630mm

VERSION	EFFECTIVE LENGTH	TOTAL LENGTH
Plain	500mm	530mm
TC	500mm	530mm
FT	500mm	530mm
Tiemann Tip	500mm	530mm

SIZES										
FG	5	6	8	10	12	14	16	18	20	22
Colour Code	Grey	Green	Blue	Black	White	Green	Orange	Red	Yellow	Violet