

Steri-Vac Set

Closed Wound Suction Unit

Closed Wound Suction Set (with Acc.-Redon Catheter PVC) | SMD 1003
Acc. Redon Drainage Catheter (PVC) | SMD 1010

- Tertutup Luka Suction Satuan
- Unidad De Aspiración Herida Cerrada
- Fechado Unidade de sucção da ferida
- Fermé Unité d'aspiration des plaies



- Used for closed wound drainage under negative pressure post-surgical patients
- Made from non-toxic, non-irritant PVC
- Radio opaque drainage catheters have multi-perforation for efficient drainage
- Option to clinicians to use one or two catheters simultaneously
- Set contains – Easy to use flexible bottle chamber with hanger, connecting tube with matching size connectors, needle with perforated catheter, one extra separate perforated catheter

CAPACITY
30ml
50ml
500ml
800ml

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 1; Master Carton of 20.

Gross: 7.2kg (14 FG) 670 x 300 x 480mm

SIZES						
FG	8	10	12	14	16	18
Colour Code	Blue	Black	White	Green	Orange	Red