

Sterimed Salem Sump Tube / Irrigation Tube

Gastro Duodenal Dual Flow Tube

SMD 600

-  Salem Sump Tube/Irrigation Tube Gastro Duodenal Dual Flow Tube
-  Salem Sumidero Tubo / Tubo De Irrigación Tubo Gastroduodenal Dual Flow
-  Salem Sump tubo / tubo de irrigação Gastro Duodenal dupla Fluxo Tubo
-  Salem puisard Tube / Irrigation Tube Gastro duodénal Dual Flow Tube



- Used to prevent the accumulation of gastric secretions and air
- Made from non-toxic medical grade PVC
- Double lumen tube throughout the length provides airvent when the stomach has been evacuated of accumulated liquid, air or small particulate matter
- Coned distal end with ten lateral eyes for faster drainage
- Universal funnel connector along with double taper connector at proximal end for easy extension
- Markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement



EFFECTIVE LENGTH	TOTAL LENGTH
119cm	122cm

Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 50; Master Carton of 600.



Gross: 24kg (16 FG)



12 to 20 FG - 680 x 480 x 520mm

SIZES					
FG	12	14	16	18	20
Colour Code	White	Green	Orange	Red	Yellow



DEHP FREE


Sterimed Levin's Tube / Stomach Tube

SMD 602

-  Levin's Tube / Stomach Tube
-  Tubo Levin's / Tubo De Estomago
-  Tubo / tubo no estômago de Levin
-  Tube / Ventre Tube de Levin



- Used for gastro-intestinal drainage
- Manufactured from non toxic, medical grade PVC
- Coned distal end is coned with four lateral eyes for efficient drainage
- Proximal end universal funnel connector
- Markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement



Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 50; Master Carton of 500 (6 to 14 FG)
- Box of 50; Master Carton of 600 (16 to 20 FG)



Gross: 16.5kg (14 FG)



6 to 14 FG - 790 x 480 x 280mm

16 to 20 FG - 680 x 480 x 520mm

EFFECTIVE LENGTH	TOTAL LENGTH
122cm	125cm

SIZES								
FG	6	8	10	12	14	16	18	20
Colour Code	Green	Blue	Black	White	Green	Orange	Red	Yellow