



Clearpore

Transparent Perforated Plastic Film Tape
SS 106 C


-  Transparan Tape berlubang Plastik Film
-  Cinta de película plástica transparente perforada
-  Fita de filme plástico transparente perfurada
-  Ruban de film plastique transparent perforé



Breathable



Water resistant



Bi-directional tearing

- Made from perforated plastic film for good breathability
- Uniformly coated with hypoallergenic adhesive
- Strong adhesion water resistant
- Painless removal
- Perfect bidirectional tearing

Packing

- Carton of 48 Boxes

SIZES & PACKING		
Size	Carton Size	Gross Weight
1/2", 1", 2" x 9.1mtr	510x465x242mm	20.50kg
3" x 9.1mtr	440x330x350mm	20.50kg
1/2", 1", 2" x 5.0mtr	400x365x240mm	12.00kg
3" x 5.0mtr	360x270x335mm	12.00kg