

Sterimed Urine Collecting Bag

-  Kantong Pengumpul Urine
-  Bolsa recolectora de orina
-  Bolsa coletora de urina
-  Sac de collecte d'urine



- Used for short and long term urine collection drained through a catheter
- 2000ml capacity bag made from medical grade soft PVC
- Non-return valve to avoid back flow of Urine
- Kink resistant 95cm tube with universal connector for easy connection to the drainage catheter.
- Graduated scale to measure urine output accurately

Product Name	Type	Product Code	Specification	Box Qty.	Master Carton	Master Carton Dimensions
Urine Bag With Hanger	With Top Outlet	SMD 504-HTP	• Specially molded handle for easy handling • Top outlet for single hand operation and efficient emptying of the bag • Capacity : 2000 ml, Tube length: 95-100cm	25 Pcs	400 Pcs	680 x 480 x 520 mm
	With Bottom (T-Type) Outlet	SMD 504-HBT	• T-type bottom outlet for single hand operation and efficient emptying of the bag			
	With Bottom Outlet	SMD 504-HBS	• Bottom outlet for single hand operation and efficient emptying of the bag			
Urine Bag Without Hanger	With Top Outlet	SMD 504-WTP	• Top outlet for single hand operation and efficient emptying of the bag • Capacity : 2000 ml, Tube length: 95-100cm	25 Pcs	400 Pcs	680 x 480 x 520 mm
	With Bottom (T-Type) Outlet	SMD 504-WBT	• T-type bottom outlet for single hand operation and efficient emptying of the bag			
	With Bottom Outlet	SMD 504-WBS	• Bottom outlet for single hand operation and efficient emptying of the bag			

Packing

- Sterile, Individually packed in Tear H.M pack.