

Sterimed Silicone Hematuria Catheter
3 WAY

Product Code for India | SMD 519 R
Product Code for Exports | SMD 519

-  Sterimed kateter hematuria silikon
-  Sterimed catéter de hematuria de silicona
-  Sterimed cateter de hematúria em silicone
-  Sterimed cathéter d'hématurie en silicone



- Specifically designed for post-TURP hemostasis and efficient hematuria management.
- High strength polymer reinforcement shaft prevents collapse under high-flow irrigation, suction, or balloon inflation pressure.
- Reinforcement ensures blockages can be cleared easily by continuous irrigation without catheter failure.
- Three-way design allows for simultaneous irrigation, drainage, and balloon inflation for optimal clot evacuation.
- Made of 100% medical-grade silicone for superior biocompatibility and patient comfort, eliminating latex allergy risks.
- Smooth silicone surface minimizes friction and tissue trauma, supporting safer and more comfortable insertions.
- X-ray opaque (radiopaque) line allows straightforward confirmation of catheter placement and intubation via imaging.
- Soft, uniformly inflated large-volume balloon sits securely against the bladder wall, preventing displacement and ensuring minimal leakage.

Packing

- Sterile, Individually packed in peelable pouch pack.
- Box of 10; Master Carton of 400 & 600

Size	Balloon Capacity
18 - 20 Fr	50 - 80 cc/mL
22 - 24 Fr	80 - 100 cc/mL



Gross: 22.0kg (600 pcs) (18FR)
Gross: 14.5kg (400 pcs) (18FR)



555 x 503 x 555mm (600 pcs)
550 x 330 x 535mm (400 pcs)