

Sterimed Urine Collecting Bag With Measured Volume Meter

Measured Volume Meter 250 ML With Bottom Outlet | SMD 510-250BS
 Measured Volume Meter 250 ML With Bottom (T-Type) Outlet | SMD 510-250BT
 Measured Volume Meter 500 ML With Bottom Outlet | SMD 510-500BS
 Measured Volume Meter 500 ML With Bottom (T-Type) Outlet | SMD 510-500BT

-  Urine Collecting Bag With Measured Volume Meter
-  Bolsa Colectora De Orina Con Medidor De Volumen
-  Urine Bag Coletando Medido com medidor de volume
-  Collecteur d'urine Sac Avec Mesuré indicateur de volume

- Used for measurement of hourly urine output
- 2000ml capacity bag made from non-toxic, non-irritant, medical grade PVC
- Rigid 250/500ml capacity bottle for accurate measurement of urine output
- Bottle directly connected to the bag to avoid contamination
- Push/Pull type bottom outlet for convenient emptying of urine bag
- Kink resistant 90cm tube with universal connector



Packing

- Sterile, Individually packed in peelable blister pack.

Product Name	Type	Product Code	Specification	Box Qty.	Master Carton	Master Carton Size
Urine Collecting Bag With Measured Volume Meter 250 ML	With Bottom (T-Type) Outlet	SMD 510-250BT	• Inlet tube length : 125 cm • Capacity of urine collection bag : 2000 ml • Bottle capacity : 250ml • With T-type Bottom outlet	5 Pcs	50 Pcs	710 x 255 x 695 mm
	With Bottom Outlet	SMD 510-250BS	• With Bottom outlet			

* Also Available in 500 ML capacity bottle.

Sterimed Leg Bag

-  Leg Bag
-  Pierna Bagsmd
-  perna Bag
-  Sac jambe



- Used for short and long term urine collection for incontinence patients
- Made from soft medical grade PVC
- Available in 750ml capacity bag, with bottom outlet for drainage, and non-return valve to prevent urine backflow.
- Graduated scale for ease of urine outflow measurement.
- Available with convenient thigh strap.

Packing

- Sterile, Individually packed in peelable blister pack.

Product Name	Type	Product Code	Specification	Box Qty.	Master Carton	Master Carton Size
Leg Bag	With Bottom (T-Type) Outlet	SMD 516-BT	• Provides freedom of mobility as the bag is attached to the leg of the patient • Capacity : 750 ml • Tube length: 30cm & 50cm • Available with leg straps • T-type bottom outlet for efficient emptying of the bag	10 Pcs	100 Pcs	570 x 340 x 360 mm
	With Bottom Outlet	SMD 516-BS	• Bottom outlet for efficient emptying of the bag			