TRANSFUSION



3 Way Stop Cock With Extension Line

SMD 707

3 Way Stop Cock With Extension Line

3 Llave De Paso Con La Línea De Extensión

3 Way Parar Galo Com a extensão da linha

3 Way robinet d'arrêt avec l'extension Ligne



· Made from non-toxic, medical grade PVC

• Has two female luer lock ports and one male luer lock port, with rotatory lock for safe connection to other devices

- Three arm handle for smooth and easy operations
- Minimum residual fluid volume
- Optional- 25/50/100/200cm lengths available in extension line

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 50; Master Carton of 500 (10cm to 50cm)
- Box of 50; Master Carton of 600 (100cm to 200cm)



Gross: 14kg (150cm)



10cm to 50cm - 790 x 480 x 280mm 100cm to 200cm - 680 x 480 x 520mm



3 Way Stop Cock

3 Way Berhenti Cock

3 Way Llave De Paso

3 Way Parar Galo

3 Way robinet d'arrêt

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 100; Master Carton of 1000.



Gross: 9kg



790 x 480 x 280mm

- · Used for pressure infusion and blood pressure monitoring
- · Made from non-toxic, medical grade PVC
- Has two female luer lock ports and one male luer lock port, with rotatory lock for safe connection to other devices
- Three arm handle for smooth and easy operations
- · Minimum residual fluid volume

www.sterimedgroup.com