

## Endotracheal Tube - Cuffed & Plain

(Nasal/Oral)

Cuffed | SMD 701 C

Plain | SMD 701 P

- Endotracheal Tube - Cuffed hidung / lisan
- Tubo Endotraqueal - Cuffed & Plain
- Tubo endotraqueal - Cuffed & Plain
- Tube endotrachéale - revers & Plaine



PHTHALATE FREE



- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade clear PVC
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Smooth inside/outside surface for ease in intubation, extubation and suctioning
- Smooth polished Murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use

### Packing

- Sterile, Individually packed in peelable pouch / blister packed.
- Box of 25; Master Carton of 200/400

Gross: 6.5kg (300 pcs)  
Gross: 12.5kg (300 pcs)

660 x 320 x 380mm (2.0 to 4.5mm)  
660 x 380 x 500mm (5.0 to 10.0mm)

#### SIZES - CUFFED

2.5 mm, 3.0 mm, 3.5 mm, 4.0 mm, 4.5 mm, 5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm, 9.0 mm, 9.5 mm, 10 mm

#### SIZES - PLAIN

2.0 mm, 2.5 mm, 3.0 mm, 3.5 mm, 4.0 mm, 4.5 mm, 5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm, 9.0 mm, 9.5 mm, 10 mm