

## Nasopharyngeal Airways

SMD 720

- Nasofaríngeo Airways
- Nasopharyngeal Airways
- nasofaríngeas Airways
- nasopharyngés Airways



PHTHALATE FREE



- Used to secure and open airways for patients in case of emergencies
- Provided with rotatable flange for easy securing and pulling out.
- Made from soft, non-irritant medical grade PVC
- Anatomically designed with smooth round edges, thin wall and kink resistance for maximum patient comfort
- E.T.O. sterile & pyrogen free. For single use only

### Packing

- Sterile individually peelable pouch packed.
- Box of 25; Master Carton of 200 pcs.

Gross: 4.2kg

470 x 240 x 350mm

### SIZES

3.0 mm, 3.5 mm, 4.0 mm, 4.5 mm, 5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm

## Guedel Airways / Oropharyngeal Airways

SMD 705

- Guedel Airways / Orofaríngeo Airways
- Cánulas De Guedel / Orofaríngeos Airways
- Guedel Airways / orofaríngeos Airways
- Guedel Airways / oropharyngés Airways



- Used to maintain an unobstructed oro-pharyngeal airway
- Made from rigid, non-toxic medical grade PVC
- Colour coded in-built bite blocks to avoid airway occlusion and biting of tongue

### Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 50; Master Carton of 400.

Gross: 8kg (Size 03)

000 to 1 - 770 x 330 x 450mm  
2 to 5 - 770 x 360 x 450mm

### SIZES

Size	Colour Code	Size	Colour Code
000	-	2	Green
00	Blue	3	Yellow
0	Black	4	Red
1	White	5	Orange