

Cytobrush /Endocervical Brush

SMD 804

-  Cytobrush / Endoserviiks Sikat
-  Citocepillo/Cepillo Endocervical
-  Cytobrush / Endocervical Escova
-  Cytobrosse / Endocervical Brush



- Used for Cervical Papanicolaou (Pap) smear screening
- Consists of Brush and the Wooden spoon Strip

Packing

- Sterile, individually packed in peelable pouch pack.
- Box/Pack of 50. Carton of 500

 Gross: 7.5kg  250 x 110 x 110mm

Infant Mucus Extractor

SMD 803

-  Infant Mucus Extractor
-  Aspirador Infantil De Mucosidades
-  Infant Muco Extractor
-  Infant mucus Extractor



- Used for aspiration of secretion from oropharynx in newly born babies
- Made from non-toxic, medical grade PVC
- Transparent chamber for visual examination of the content
- Round open tip for atraumatic intubation, and proximal end with funnel connector for easy connection to the suction equipment
- Optional: available with filter and without filter models

Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 25; Master Carton of 400.

 Gross: 13.5kg  670 x 300 x 480mm

SIZES	
Chamber Capacity	25ml
Catheter Size	10 FG