

Oxygen Mask

SMD 704

- Oxygen Mask
- Máscara De Oxígeno
- Máscara de oxigênio
- Masque d'oxygène



- Used for efficient oxygen administration to patient
- Made from soft, non-toxic, medical grade PVC
- Elastic strap and aluminum nose clip for proper adjustment of mask
- Multi-channel tube ensures regular supply of oxygen even if the tube is accidentally kinked
- Proximal end of tube provided with universal connector for easy connection with oxygen source

Packing

- Non- Sterile, individually packed in a polybag.
- Box of 10; Master Carton of 200.

Gross: 24.4kg (Adult) 680 x 480 x 520mm

SIZES	TUBE LENGTH
Adult	210cm (7ft)
Child	210cm (7ft)

Nasal Oxygen Cannula

Twin Bore Oxygen Delivery Set

SMD 703 A

- Nasal Oxygen Cannula Twin Bore Oxygen Delivery Set
- Canula Nasal De Oxigeno Entrega Con Doble Juego
- Nasal cânula de oxigênio Bore Twin Set entrega de oxigênio
- Canule nasale d'oxygène Bore Twin Set oxygène de livraison



- Used for efficient oxygen administration to patients
- Made from soft, non-toxic, medical grade PVC
- Twin nasal tips ensure equal volume of oxygen to both nasal passages for greater patient comfort
- Multi-channel tube ensures regular supply of oxygen even if the tube is accidentally kinked
- Proximal end of tube provided with universal connector for easy connection with oxygen source

Packing

- Non-Sterile, individually packed in a polybag.
- Box of 20, Master Carton of 200.

Gross: 13.8kg 670 x 300 x 480mm

SIZES	TUBE LENGTH
Adult	210cm (7ft)
Child	210cm (7ft)
Infant	210cm (7ft)