TAPES & DRESSING



Steripore / Steripore+ Microporous Tape

SS 106

Tape Mikro

Cinta Microporosa

Tape microporosa

bande microporeux





Soft and Hypoallergenic



Easy to cut and tear



Suits contoured surfaces



Highly porous – allowing the skin to breathe



Firm yet easy to remove



Suitable for sensitive skin

- Used for fixation of dressing, tubing, osteomy pouches etc on patients
- Made from micro-porous non-woven fabric
- Uniformly coated with hypoallergenic, non-toxic, non-irritant adhesive
- Easy to peel off without causing discomfort or skin irritation
- · High porosity allows skin to breathe, and high resistance to water and vapour prevents the tape from peeling off
- Suitable on sensitive body parts, for infants and adults
- Manufactured as per British Pharmacopia (BP) guidelines

Packing

· Carton of 60 Boxes

SIZES & PACKING		
Size	Carton Size	Gross Weight
1/2", 1", 2" x 9.1mtr	440x405x285mm	15kg
3" x 9.1mtr	510x305x335mm	15kg
1/2", 1", 2" x 5.0mtr	390x340x285mm	10.1kg
3" x 5.0mtr	430x270x335mm	10.1kg

www.sterimedgroup.com 46