



EN ISO 13485:2016

STERIMED P.P.E. KIT

Personal Protective Equipment Kit

Specially Developed to Protect against COVID-19 Infection



Coveall Suit

1. DRDE approved Coverall suits with hood and long shoe cover made of DRDE approved SMS nonwoven fabric laminated with GSM 75/80 with seam sealing tape.

**Nitrile gloves
Megacare brand
Gloves cobalt color**

2. Disposable Cobalt Color Nitrile Gloves with texture fingers high standard of Megacare brand Approved EN ISO 13485:2016



Face Shield

3. Sterimed face shield of plastic professional face with anti-fog, anti-glare performance, heavy duty and most durable with the adjusting band according to different sizes.

OR



We like to introduce ourselves as DRDE approved supplier of PPE kits coverall. Above & Below is the brief of components of PPE kits.



EN ISO 13485:2016

STERIMED P.P.E. KIT

Personal Protective Equipment Kit

Specially Developed to Protect against COVID-19 Infection



Plastic Goggles

Branded Goggles heavy and durable frame of Goggles with good seal with the skin of face, covering the eye and the surrounding areas and even accommodating for prescription glasses with the adjusting band according to different sizes.

4. Mask - a triple layer (3-ply) medical disposable Mask fluid-resistant for droplets protection with machine sealed and stitched with elastic band to fit in the ear along with the nose pin to seal nose bone.



**Disposable Mask
(3 Layers)**



Shoe Cover

5. Made of poly propylene nonwoven with lamination and open side sealed/ stitched with elastic band, Which prevents any kinds of penetration.

6. The Waste / Disposable Bag are used to Packed, when remove Personal Protective Equipment.



Waste/ Disposable Bag

Manufactured by :

Sterimed Medical Devices Pvt. Ltd.

Plot No. 211A, Sector-16, HSIIDC, Bahadurgarh, Jhajjar-124507, Haryana (INDIA)

Customer care No.: +91-11-48880000

E-mail : info@sterimedgroup.com Website : sterimedgroup.com.

Mfg. Lic. No. MFG/MD/2018/000086

