



P.E. Pressure Monitoring Line

Spiral

SMD 715

-  P.E. Pressure Monitoring Line Spiral
-  P.E. Línea Espiral De Monitoreo De Presión
-  P.E. Linha de Monitorização da Pressão
-  P.E. Surveillance de la pression de ligne

- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from Kink resistant PE Tubing
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Internal Diameter 1.0 mm, Outer Diameter 2.0 mm

Packing


- Sterile, Individually packed in peelable pouch & blister pack.
- Box of 50; Master Carton of 500



Gross: 13.5kg



790 x 480 x 280mm



SIZES			
Length	110cm	150cm	200cm
Prime Volume	0.9ml	1.40ml	1.80ml

Pressure Monitoring Line

SMD 706

-  Pressure Monitoring Line
-  Línea De Monitoreo De La Presión
-  Linha de Monitorização da Pressão
-  Surveillance de la pression de ligne



- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from non-toxic medical grade PVC
- Can withstand a pressure of 580 PSI
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Optional lengths available

Packing

- Sterile, Individually packed in peelable pouch & blister pack.
- Box of 50; Master Carton of 500



Gross: 14.5kg



790 x 480 x 280mm

SIZES			
Length	110cm	150cm	200cm
Prime Volume	0.9ml	1.40ml	1.80ml