

## Stericiser

Breathing Exerciser  
SS 1031

- Berolahraga Pernapasan
- Ejercitador Respiratorio
- respiração Exerciser
- respiration exerciseur

- Used as a breathing exerciser to increase lung capacity in patients with decreased lung function
- Made from rigid, non-toxic medical grade PVC
- Specially designed with 3 stage chambers to provide more efficient exercise to patient at a step up basis
- Device comprises a base and a central part with a 3 stage chamber, each containing spheres of different sizes, weight and colour
- Provides lung exercise on a step-up basis, as the spheres rise proportionately to the breathing effort exerted by the patient
- Can be dis-assembled into parts for cleaning and disinfection
- Non-sterile

### Packing

- Sterile / Disposable / Individual Printed Box Packing.
- Box of 1/Master Box of 50.

Gross: 11.3kg 710 x 5000 x 350mm



## Steri-Vac Set

Closed Wound Suction Unit  
SMD 1003

- Tertutup Luka Suction Satuan
- Unidad De Aspiración Herida Cerrada
- Fechado Unidade de sucção da ferida
- Fermé Unité d'aspiration des plaies



PHthalate FREE

- Used for closed wound drainage under negative pressure post-surgical patients
- Made from non-toxic, non-irritant PVC
- Radio opaque drainage catheters have multi-perforation for efficient drainage
- Option to clinicians to use one or two catheters simultaneously
- Set contains – Easy to use flexible bottle chamber with hanger, connecting tube with matching size connectors, needle with perforated catheter, one extra separate perforated catheter

### Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 1; Master Carton of 20.

Gross: 7.2kg (14 FG) 670 x 300 x 480mm

SIZES						CAPACITY
FG	10	12	14	16	18	500 ml
Colour Code	Black	White	Green	Orange	Red	600 ml