

## Yankaur Suction Set

With Standard Handle

SMD 712

-  Yankaur Suction Set With Standard Handle
-  Juego De Aspiración Yankaur Con Asa Standard
-  Yankaur Conjunto de Sucção Com Punho Padrão
-  Yankaur aspiration avec poignée standard

- Used for removal of secretion and blood pass during surgical procedures
- Made from non-toxic, non-irritant PVC
- Flexible kink resistant tube provides uninterrupted flow during suction
- Soft connectors at both ends of the tube for easy connection
- Available in Standard / Round Tip Handle
- Suitable for use with or without vent control facility



### Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 10; Master Carton of 100

 Gross: 18kg

 740 x 425 x 425mm

SIZES	
Tube Length	2500mm (2.5mtr.)
Tube Diameter	9mm

## Surgical Blade With Handle

SMD 1006

-  Pisau bedah Dengan Handle
-  Hoja Quirúrgica Con Asa
-  Lâmina cirúrgica Com Handle
-  Lame chirurgicale Avec Poignée

- Disposable, Sterile, ready to use
- Fitted on handle No. 3 or 4



### Packing

- Pack of 20. Carton of 200pcs.

 Gross: 1.78kg

 152 x 229 x 254mm

### SIZES

10, 11, 12, 15, 20, 21, 22, 23, 24