

## Surgeon Cap

SMD 1204

## Nurses Cap

SMD 1206

- Gorro De Enfermeras
- Gorro De Cirujano
- Cap cirurgiã / tampão das enfermeiras
- Chirurgien CAP / infirmières Cap

- Used to protect against contamination from the head/hair
- Made from non-woven fabric
- Provided with stitched interlocking and elastic strip on the back for proper grip
- Optional: Available in Blue /Green colour

### Packing

- Individually packed.
- Box of 100; Carton of 1800pcs.

Gross: 9.51kg 660 x 432 x 508mm

## Surgical Gown

SMD 1202

- Bedah Gown
- Bata Quirúrgica
- Vestido cirúrgico
- casaque chirurgicale

- Used as a protective barrier to prevent cross-contamination
- Made from soft non-woven fabric with stitched interlocking
- Light weight, with soft elastic cuff and tie straps for better fitting and complete body wrap

### Packing

- Individually packed in polybag
- Box of 50.
- Master Box of 100kg.

Gross: 15kg 533 x 533 x 558mm

## Shoe Cover

SMD 1210

- Shoe Cover
- Cubierta De Zapatos
- Shoe Cover
- Couvre-chaussure

- Used to protect against contamination from footwear
- Made from non-woven fabric/polyethylene, on a fully automatic ultrasonic sealing machine
- Elastic strip for proper grip

## Surgeon Face Mask

SMD 1208

- Dokter bedah Masker Wajah
- Mascara Facial De Cirujano
- Cirurgião Máscara Facial
- Visage Chirurgien Masque

- Used as a protective barrier to prevent cross-contamination
- Made from non-woven fabric and available in two-layers and three - layers configuration
- Ultrasonically welded layers for efficient bacterial filtrations
- Provided with a nose clip to prevent fogging of spectacles and to maximize filtration efficiency

### Packing

- Individually packed.
- Box of 100; Carton of 4000.

Gross: 19.57kg 660 x 432 x 508mm

## Surgical Gloves

SMD 1000

- Sarung tangan bedah
- Guantes Quirurgicos
- luvas para cirurgia
- Gants chirurgicaux

- Used for protection during general medical /surgical procedures
- Made from natural latex rubber
- Firm beaded cuffs, micro rough finish for better grip
- Bent fingers for greater comfort
- Lightly powdered

### Packing

- Carton of 1000 pairs. Sterile

Gross: 28.3kg 483 x 660 x 457mm

### SIZES

6.0, 6.5, 7.0, 7.5, 8, 8.5

### Packing

- Individually packed in polybag
- Carton of 2000 pairs. Pack of 25 pairs.

Gross: 15.4kg 737 x 584 x 432mm

