





Intubating Stylet

SS 753

-  intubasi stylet
-  intubación estilete
-  intubating Stylet
-  stylet d'intubation

- Used along with Endotracheal Tube to aid during intubation
- Smooth plastic coated with flexible malleable aluminum design, allows the stylet to be easily adapted to any desired shape
- Soft distal tip to reduce the possibility of intubation damage
- Smooth surface for ease in introduction and removal from ET tube
- Clean packed

Packing


- Individually pouch packed.
- Inner Box of 10 pcs.
- Master Carton of 200 pcs.



Gross: 6.6kg



420 x 325 x 280mm








SIZES

6 Fr/Ch (2.0 mm)
10 Fr/Ch (3.3 mm)
14 Fr/Ch (4.7 mm)

Visionary Bougie Vented

SS 755

-  Visioner Bougie Vented
-  Visionario Bougie ventilado
-  Visionário Bougie Vented
-  Visionnaire Bougie ventilé



- Provides rigidity to endotracheal tubes to facilitate endotracheal intubation
- Ultra smooth surface for low friction
- Atraumatic angled tip for ease of intubation and to minimise risk of patient trauma
- Marking provided along the length to facilitate identification of tube position
- Angled tip for ease of intubation
- Vented lumen for oxygen feed

Packing

- Sterile, individually packaged in easy to open sachets
- Master Carton of 10 pcs.



Gross: 0.5kg



800 x 75 x 60mm

SIZES

FR	6	10	15
----	---	----	----



Atraumatic Smooth Tip



Removable Soft Connector