





Infusion Set

SMD 1100

-  Perangkat Infus
-  Set de infusion
-  Conjunto de infusão
-  Set de perfusion



- Sharp piercing spike for easy insertion in I.V Container
- Cylindrical collapsible drip chamber to visualize the flow rate
- Disc type fluid filter to filter any particulate matter in the I.V fluid
- Super smooth kink resistant tubing with efficient roller controller for accurate and unrestricted flow
- Flash ball type injection port for extra medication
- Infusion set with drop rate of 20 drops/ml

Packing

- Sterile, individually packed in poly pouch.
- Special puncture resistant packing
- Inner poly bag 25 pcs.
- Master shipper 600 pcs.



Special puncture resistant packing



Gross: 16.9kg



635 x 485 x 380mm